




# Action Groups Guide

The Action Groups feature of AirView filters your patients into actionable groups based on certain criteria. It helps you manage your patients efficiently by focusing your attention specifically on patients that need your assistance.

Action Groups allow you to manage your patient population efficiently. For selected organisations, Action Groups will be activated in your AirView account. Patients are added/removed from a group based on having met a certain criteria. Additionally, you can exempt the patient from a certain parameter.

To navigate, from the Patients menu, select Action Groups.




#### Therapy issues Usage issues No data received Ongoing therapy follow up Ongoing usage follow up Identifies patients who have experienced Identifies patients who have demonstrated Identifies patients whose data have not been Identifies patients who are experiencing Identifies patients who have demonstrated a higher than expected mask leak, AHI or CAI a trend in low usage on their device in their transmitted for 3 days in their first 30 days higher than expected mask leak, AHI or CAI trend in low usage on their device after their in their first 30 days of therapy. first 30 days of therapy. of therapy and for 10 days after the first 30 after their first 30 days of therapy. first 30 days of therapy. days.



### Action Groups: auto-add rules


AirView automatically adds patients into Action Groups when they meet the relevant criteria and removes them when they no longer meet those criteria. Patients can also be manually removed from Action Groups. Please see pages 5 and 6 for details on auto and manual removal.



- Leak is always 24 L/min and is not changed by preferences: Patient Leak Threshold
- No data triggers inclusion in the No data received group, but does not add the patient to other groups
- To disable the auto-add function for individual patients, go to Patient profile > Preferences > Edit.



# Features available inside each group





## Features available inside each group


#### View recent logs

From the Action menu, click View recent logs to display the reason why the patient was added into the Action Group. You can find all Action Group logs when you click View all Action Group logs.



### View reviewed patients

Click View reviewed patients to display all patients that have been reviewed in the last 48 hours.



#### Mark as reviewed

Once you have actioned any patient issues (e.g. called the patient, adjusted settings), you can remove the patient from the Action Group.

From the Patients > Action


From the Patients > Action Groups menu, select the patient and click Mark as Reviewed.

You can also select multiple patients or select all patients in a group.



### Restore reviewed patients


Select the patient and click **Restore selected to group**, to restore a patient back to the Action Group.





### Action Groups: auto-remove rules


Auto-remove rules simplify the management of Action Groups. Patients are automatically removed from an Action Group when they no longer meet any of the criteria for inclusion. This change of status is recorded in the patient groups log. Patients who are auto-removed do not appear in the **Mark as reviewed** list. Patients who are manually removed from a group appear in the **Reviewed patients** list.





## Excluding a patient from Action Groups

**Example:** Joe has been auto-enrolled in the **Therapy issues** Action Group due to high leak. I know he is doing well with therapy in spite of his leak issues so I will remove him from the group by following steps 1 to 4.



- The manual removal option is only available to organisations with Action Groups.
- If a patient was previously excluded from a group for a specific metric (e.g. leak, usage, etc), upon inclusion they can be triggered the next day.