# NV Elbow

F20

Non-vented elbow

ar	كوع غير منفس
bg	Невентилирано коляно
cz	Neventilované kolínko
da	Ikke-ventileret vinkelstykke
de	Kniestück ohne Luftauslassöffnungen
el	Γωνιακό εξάρτημα χωρίς εξαερισμό
en	Non-vented elbow
es	Codo sin orificios de ventilación
et	Mitteventileeritav torupõlv
fa	زانوی ماسک بدون تهویه
fi	Ilma-aukoton kulmakappale
fr	Coude non ventilé (sans fuite intentionnelle)
hr	Neventilirano koljeno
hu	Nem szellőző könyök
is	Olnbogi án loftopa
it	Gomito non-vented
lt	Neventiluojiama alkūnė
lv	Izliekums bez ventilācijas atverēm
nl	Niet-geventileerde elleboog
nn	Ikke-ventilert kneledd
pl	Kolanko bez wylotów powietrza
pt	Cotovelo não ventilado
ro	Cot neventilat
ru	Невентилируемый угловой патрубок
sl	Neventilirano koleno
sv	Oventilerat knärör
tr	Ventilsiz dirsek

## User Guide

English | Français canadien | Español de América | Português do Brasil | Français | Español  
Português | Български | Hrvatski | Český | Dansk | Nederlands | Eesti | Suomi | Deutsch  
Ελληνικά | Magyar | Íslenska | Italiano | Latviski | Lietuvių k. | Norsk | Polski | Română | Русский  
Slovenščina | Svenska | Türkçe | فارسی | العربية

**ResMed Pty Ltd**  
1 Elizabeth Macarthur Drive  
Bella Vista NSW 2153 Australia

See ResMed.com for other ResMed locations worldwide. AirFit is a trademark and/or registered trademark of the ResMed family of companies. For patent and other intellectual property information, see ResMed.com/ip. © 2020 ResMed. 638245/2 2020-11

ResMed.com

Rx Only

CE 0123

638245

NV Elbow F20  
Non-vented elbow

## ENGLISH

The NV (Non-Vented) Elbow F20 replaces the existing vented elbow in the AirFit™ F20 mask and converts it to a non-vented mask system.

### Intended use

The NV (Non-Vented) Elbow F20 converts the AirFit F20 into a non-invasive mask used for channeling airflow (with or without supplemental oxygen), intended to be used with active-exhaust-valve ventilators, to provide ventilatory assistance to patients with respiratory insufficiency and respiratory failure.

The NV Elbow F20 (when used with the AirFit F20 mask) is:

- to be used by patients weighing >66 lb (30 kg) requiring non-life support ventilatory assistance.
- intended for single-patient reuse in the home and hospital/institutional environment.

### GENERAL WARNINGS

Refer to the AirFit F20 User Guide for additional warnings related to using the mask. The following warnings relate to the use of the NV Elbow F20.

- The NV Elbow F20 does not contain an exhaust vent or safety valve. The non-vented elbow must be used with breathing circuits or PPV devices which provide their own method of venting expired or supplemental gases.
- The mask must be used under qualified supervision for users who are unable to remove the mask by themselves or when the patient is at risk of not responding to severe hypoxia or hypercapnia by arousal. The mask may not be suitable for those predisposed to aspiration.
- The mask should only be used with ventilators recommended by a physician or respiratory therapist.
- The mask should not be worn unless the ventilator system is turned on and operating properly.
- The technical specifications of the mask are provided for your clinician to check that they are compatible with the ventilator. If used outside specification or if used with incompatible devices, the seal and comfort of the mask may not be effective, optimum therapy may not be achieved, and leak, or variation in the rate of leak, may affect the ventilator function.
- Visible criteria for product inspection: If any visible deterioration of a system component is apparent (cracking, discoloration, tears, etc.), the component should be discarded and replaced.

**Note:** Any serious incidents that occur in relation to this device should be reported to ResMed and to the competent authority in your country


### Using the elbow

ResMed has tinted the NV Elbow F20 blue to identify it as a non-vented accessory. The elbow has a standard female 22mm tapered fitting for external connection to the ventilator. Refer to the Technical specifications section in this user guide for information about mask leak, dead space and resistance.

### Assembly

1. If the vented elbow is still assembled to the mask, squeeze the side buttons and pull the elbow away from the frame.

2. Attach the NV Elbow F20 to the mask by pushing the elbow into the frame, ensuring it clicks in and is firmly secured.

## Disassembly

Remove the NV Elbow F20 from the mask by squeezing the top and bottom side buttons and pulling away from the frame.

## Cleaning

**Daily/After each use:** Refer to the AirFit F20 User Guide. The cleaning instructions for the elbow apply to the NV Elbow F20. Discard if heavily soiled.

## Reprocessing between patients

The NV Elbow F20 is intended for single-patient reuse only. This component cannot be reprocessed and must be replaced between patients. The AirFit F20 mask can be reprocessed, refer to the AirFit F20 reprocessing instructions.

## Technical specifications

The following specifications apply to the NV Elbow F20 with the AirFit F20 mask.

**Compatible devices:** Astral

**Deadspace:** Physical dead space (mask and elbow) using the large cushion is 248 mL.

**Therapy pressure:** 3 to 40 cm H<sub>2</sub>O

**Pressure-flow**

3 – 20 cm H<sub>2</sub>O: < 6 L/min

20 – 40 cm H<sub>2</sub>O: < 12 L/min

**Resistance:** Drop in pressure measured (nominal)

at 50 L/min: 0.2 cm H<sub>2</sub>O

at 100 L/min: 0.6 cm H<sub>2</sub>O

**Materials:** Polycarbonate

**Environment conditions, mask gross dimensions and service life:** Refer to the AirFit F20 User Guide

## Storage and disposal

Refer to the AirFit F20 User Guide for storage and disposal.

## Symbols

The following symbols may appear on your product or packaging:

**NV** Non-vented mask   **MD** Medical device   **⚠** Indicates a Warning or Caution and alerts you to a possible injury or explains special measures for the safe and effective use of the device   **R<sub>x</sub> Only** Prescription only (In the US, Federal law restricts these devices to sale by or on the order of a physician.)

See symbols glossary at [ResMed.com/symbols](https://www.resmed.com/symbols).

## Warranty

Refer to the AirFit F20 User Guide for warranty information.

## Further information

For the most up to date information and instructions in other languages, go to [ResMed.com/downloads/masks](https://www.resmed.com/downloads/masks).

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