


ResMed


The control panel and user interface



 Inova Labs
3500 Comsouth Drive,
Suite 100
Austin, TX 78744 USA


Distributed by: See ResMed.com for ResMed locations worldwide. ActiveVox, LifeChoice and PULSE-WAVE are trademarks and/or registered trademarks of the ResMed family of companies. For patent information, see ResMed.com/ip. Product availability may vary across regions. Specifications may change without notice. © 2017 ResMed Ltd. 1019963/1 2017-06

ResMed.com



ResMed

LifeChoice **Activox**[™] 4L



Quick set up guide

How to charge your POC


Please ensure you fully charge the internal battery before using your device for the first time. If you purchased an external battery, it should also be fully charged before use.



1 To charge the internal battery, connect the AC power supply to the device and plug it into an electrical socket. A plug icon  will appear on the display screen to show that your portable oxygen concentrator (POC) is connected to an external power source.



2 To charge the external battery, first connect the AC power supply to an electrical socket (see step 1). Then connect the external battery to the external battery port on the device (the flat side of the external battery cable should be facing upwards).



3 A 'Charging Internal' message and the 'INT' battery bar on the display screen show that the internal battery is being charged. The same indicators with 'EXT' mean the external battery is being charged. When the battery is fully charged, the message will disappear from the display and the battery bar will indicate the charge remaining.

Note: It can take **up to 4 hours** to charge a depleted internal battery and **up to 2 hours** to charge a depleted external battery. By design, the system will always charge the internal battery first and then charge the external battery. The system will deplete the external battery first to preserve the internal battery for as long as possible.

Charging tips

- We recommend keeping your device plugged into the mains power supply even when it is not in use and charging the external battery as often as possible. This helps extend the overall life of the rechargeable batteries and only consumes a very small amount of electricity. The batteries will not over-charge.
- Please ensure you test your vehicle's charging capabilities before taking a trip. Some vehicles may not have sufficient power output to the DC power supply* to simultaneously operate and charge the device and/or external battery.
- The fans on your POC will turn on to remove heat that is generated during operation and battery charging.

| Setting | 1 | 2 | 3 | 4 |
|-----------------------------------|------------|-----------|-----------|-----------|
| Internal battery runtime** | 10hr 15min | 8hr 15min | 5hr | 4hr |
| External battery runtime** | 4hr 45min | 3hr 45min | 2hr 45min | 2hr 25min |
| Total runtime** | 15hr | 12hr | 7hr 45min | 6hr 25min |

*For in-vehicle charging on the go **Hours are approximate and based on breaths per minute.

How to use your POC



- 1 Turn on your POC by briefly pressing and releasing the power button. The display screen and green power indicator will light up.

Important: do not put the cannula in your nose until the unit is turned on.



- 2 Press the + or – button to select the pulse setting prescribed by your doctor. The setting number (1-4) will appear on the bottom left-hand corner of the display screen, under 'Setting'.



- 3 Plug in the nasal cannula then put the cannula prongs into your nostrils and breathe normally. When you take a breath, the device will deliver a gentle dose of oxygen and the PULSE-WAVE™ symbol  will flash on and off.



- 4 Turn off your POC by pressing and holding the power button until you hear a beep.

Note: You might trigger the “no breath” alarm if you remove the nasal cannula while the device is on or breathe through your mouth for an extended period of time. To turn off the alarm, simply reconnect the cannula or breathe through your nose.

Device use and storage

- Your POC must only be used in well ventilated locations. Make sure that the inlets and outlets are not obstructed (e.g. by clothing...).
- Your POC is designed to be operated and stored in its carrying case at all times.
- Your POC and batteries may be damaged if they are exposed to extremes of temperature (e.g. in a parked vehicle). You should also keep your device away from environments with high levels of moisture, dust or sand.