



AirFit™ F30i
Full face mask



This is for the active sleepers,
who go to bed when it suits them.
For the night-time nomads
who twist and turn
then spring up to take on the day.



Tube-up freedom in a full-face mask


Introducing the AirFit F30i, the tube-up, full-face mask that allows patients to twist and turn in bed, sleep in any position¹, and read and relax with ease. The ultra-compact cushion is curved to fit comfortably under the nose and creates a stable seal.¹

AirFit F30i combines sensations of visual and physical freedom with real freedom of movement.




A gentle, low-key approach

The AirFit F30i is designed to overcome common barriers to therapy with a full-face mask. Claustrophobia, inconvenience, personal appearance, discomfort²: the F30i mask is engineered to tackle these challenges head on.





Product codes



AirFit F30i - mask system			
Size	Small cushion / Small conduit frame	Medium cushion / Standard conduit frame	Wide cushion / Standard conduit frame
EUROPE*	63310	63312	63313

***Europe:** English, German, French, Italian, Dutch, Swedish, Norwegian, Finnish, Danish, Spanish, Portuguese, Czech, Polish and Arabic.

A welcoming start

The AirFit F30i mask creates good conditions for new users to get used to sleep apnoea therapy. The design focus on physical, visual and psychological freedom helps to support natural sleep and facilitates acceptance.


Sizing tips

AirFit F30i has four interchangeable cushion sizes

- small, medium, small wide and wide. The conduit frame is available in three different sizes: small, standard and large.

The modular design of the AirFit F30i means each size cushion easily clicks into the different mask conduit frames, so you can swap between them with ease for more flexibility.

We recommend that you refer to the mask fitting template to gauge which size is the most suitable for each of your patients.



Patient fitting guide for AirFit F30i



- With one hand holding the bottom of the cushion and the other holding the conduit frame, position the mask in front of the face.



- Place the cushion under the nose and against the face. With the ResMed logo and headgear facing up, pull the headgear and conduit frame over the head.



- Bring the lower headgear straps under the ears and attach the magnetic clips³ to the frame connectors.



- Undo the fastening tabs on the upper headgear straps and pull evenly.



- Repeat with the lower headgear straps.
- Ensure the patient achieves a comfortable fit.



- Squeeze the side buttons on the elbow and detach from the frame. Connect the air tubing from the device to the elbow, then reattach it to the frame.

Refer to the user guide for further instructions.



Additional options

- For people with long hair, we recommend feeding the hair through the back of the headgear, as pictured. If tied in a ponytail, it should sit high enough to fit through the back of the headgear.
- Refer to the user guide for further instructions and complete labeling information including indications, contraindications, warnings and precautions.



This is for them.

1 ResMed external clinical study of 51 current full face CPAP mask users, conducted between 24/10/2018 - 20/11/2018, 2018 Data on file; ID A4484728.

2 Wimms A. et al., Impact of a New Nasal Pillows Mask on Patients' Acceptance, Compliance, and Willingness to Remain on CPAP Therapy. Sleep Disord, 2016. 2016: p. 6713236.

3 Use of masks with magnetic components is contraindicated in patients with the following pre-existing conditions: a metallic haemostatic clip implanted in your head to repair an aneurysm and/or metallic splinters in one or both eyes.