

AirCurve™ 11 VAuto Prescription Form

Indication for use:

The AirCurve 11 VAuto device is indicated for the treatment of Obstructive Sleep Apnoea (OSA) in patients weighing more than 66 lbs / 30 kg and more than 30 lbs / 13 kg in CPAP and S modes. The AirCurve 11 VAuto device is intended for use in the hospital and home.



Patient details

Name	Date of birth (DD-Mon-YYYY)	Phone
Address		
Medicare/NHI number	IRN	DVA patient File number

Therapy mode

Therapy Mode	VAuto	S	CPAP
--------------	-------	---	------

VAuto and S therapy settings

Parameter	Prescribed value				Unit	Device default	Device range
VAuto mode							
Min. EPAP:					cmH ₂ O	4	4 to 25, in increments of 0.2
Max. IPAP:					cmH ₂ O	25	4 to 25, in increments of 0.2
Pressure Support:					cmH ₂ O	4	0 to 10, in increments of 0.2
S mode							
IPAP:					cmH ₂ O	10	4 to 25, in increments of 0.2
EPAP:					cmH ₂ O	4	2 to "IPAP", in increments of 0.2
Backup Rate:	On (10BPM)	Off				Off	
Easy-Breathe:	On	Off				Off	
Rise-Time:					ms	300	Min / 150 to 900, in increments of 50
Ti Range, Trigger and Cycle (VAuto and S mode)							
Ti Max:					sec	2	0.3 to 4, in increments of 0.1
Ti Min:					sec	0.3	0.1 to "Ti Max", in increments of 0.1
Trigger:	Very Low	Low	Medium	High	Very High	Medium	
Cycle:	Very Low	Low	Medium	High	Very High	Medium	

CPAP therapy settings

Parameter	Prescribed value	Unit	Device default	Device range
Set pressure:		cmH ₂ O	10	4 to 20, in increments of 0.2

OFFICE USE	
Check 1	Check 2

AirCurve 11 VAuto Prescription Form

Supplemental oxygen settings

Does the patient need supplemental oxygen with their device? **Yes** **No**

Parameter	Prescribed value	Unit	Device default	Device range
Flow rate:		L/min	–	Available only when supplemental oxygen selected is "Yes". Up to 15L/min in all modes

Comfort settings

Ramp:	On	Off			
Parameter	Prescribed value			Device default	Device range
Ramp time:				20	Available only when ramp time selected is "On". 5 to 45, in increments of 5
Start EPAP:				4	Available in VAuto and S mode only. "Min device pressure for the therapy mode" to "EPAP", in increments of 0.2
Start pressure:				4	Available in CPAP mode only. 4 to "Set pressure", in increments of 0.2
EPR:	On (full time)	On (Ramp only)	Off	On (Ramp only)	Available only in CPAP mode
EPR level:	1	2	3	1	Available only in CPAP mode

OFFICE USE	
Check 1	Check 2

Mask type:

Full face Nasal Pillow Other

If other, please specify

Prescriber Details

Name Phone

Address Email

Provider number

Signature Date (DD–Mon–YYYY)