



**ResMed**

## Homecare Stand

The homecare stand is a portable, space-saving solution for the hospital or home. It protects against everyday hazards with added mobility, security and convenience.





## Space-saving convenience

The ResMed homecare stand is an ideal solution where space is limited and mobility is a must.

Its minimalist design offers maximum stability when stationary. Its light weight and built-in handle makes it easy to move from room to room. An integrated storage basket keeps patient accessories and supplies tidy and at hand.



Homecare stand shown with Astral™

## Key features

- Mobile stand with a built-in handle and 360° swivel wheels make moving between different environments a hassle-free experience.
- Secure frame and top mount enable you to attach a device to the stand, protecting it against day-to-day knocks, spills and accidental loss.
- Convenient basket, central stand and cable hooks reduce clutter. Accessories and components remain within easy reach.

## Technical specifications

<b>Height</b>	100 cm (39.4")
<b>Base diameter</b>	55.9 cm x 55.9 cm (22" x 22")
<b>Shipping weight</b>	17 kg (37.5 lb) approx.
<b>Load limit</b>	56.7 kg (125 lb)

## Product codes

Homecare stand without interface	<b>27963</b>
Astral top plate including stand interface	<b>27956</b>
IV pole kit (includes IV pole and crossbar to attach to the stand)	<b>27948</b>
Circuit support arm (requires IV pole kit to attach to the stand)	<b>24917</b>

## Spare parts

Astral top plate spare part (without stand interface)	<b>27961</b>
Homecare stand screws spare kit	<b>27960</b>
Homecare stand casters spare part	<b>27962</b>