

# Software release notes

English

## AirMini™ by Resmed App 1.9

### Product description

AirMini by Resmed is a mobile application used by patients and clinicians. Patients can operate their prescribed Resmed AirMini device, view collected data and change comfort settings. Clinicians can also operate the device and change the patient's prescription settings.

### What's new

- Bug fixes and performance enhancements.

### Deferred issues

iOS™ and Android™

- **Problem:** In some cases, dates and time in the AirMini app, such as therapy session start and end time, are not in the local time zone.  
**Solution:** There is currently no solution to this problem.
- **Problem:** In rare cases, AirMini data may not upload to AirView™.  
**Solution:** If you, or your healthcare provider, are missing any AirMini data, contact Resmed support.

## Supported operating systems

- Android operating system version 9 and newer
- iOS 10.2 and newer

## Compatible machines

- The AirMini by Resmed app is compatible with the AirMini device.

## Supported languages

- |                   |                           |
|-------------------|---------------------------|
| • Bulgarian       | • Japanese                |
| • Croatian        | • Korean                  |
| • Czech           | • Norwegian               |
| • Danish          | • Polish                  |
| • Dutch           | • Portuguese (Brazil)     |
| • UK English      | • Portuguese (Portugal)   |
| • US English      | • Romanian                |
| • Finnish         | • Slovakian               |
| • French (France) | • Slovenian               |
| • German          | • Spanish (Latin America) |
| • Greek           | • Spanish (Spain)         |
| • Hungarian       | • Swedish                 |
| • Indonesian      | • Turkish                 |
| • Italian         |                           |

## Still having problems?

Contact support:

AirMiniAppSupport@resmed.com



ResMed Ltd. 1 Elizabeth Macarthur Drive Bella Vista NSW 2153 Australia

See Resmed.com for other Resmed locations worldwide. Android is a trademark of Google Inc. iOS is a trademark or registered trademark of Cisco in the U.S. and other countries and is used under license. Bluetooth is a registered trademark of Bluetooth SIG, Inc. AirMini and AirView are trademarks and/or registered trademarks of the Resmed family of companies. For patent and other intellectual information, see [www.resmed.com/ip](http://www.resmed.com/ip).

© Resmed Ltd. 108373/1 2026-7

Resmed.com