

**ResMed**

Actual Size: 140 mm (Width)  
x 175 mm (Height)

Ruler is for scaling purposes.

How to use the fitting template:


1. Please ensure the template is printed to scale
2. Cut out along dotted line
3. Refer to diagram on template to determine appropriate cushion size



Nasal mask



Nasal mask



Sample size check

10 mm 20 mm 30 mm 40 mm 50 mm 60 mm 70 mm 80 mm 90 mm 100 mm 110 mm 120 mm 130 mm 140 mm 150 mm 160 mm 170 mm