




Mirage Quattro™

FULL FACE MASK

Quick Fitting Guide

Svenska




Components Card

Mirage Quattro Full Face Mask

XS* 61208 S* 61209 M* 61210 L* 61211

Elbow assembly (elbow, valve, clip, swivel) / Bágmontering (báge, klaff, klaffclips, svängtapp) / Vinkelsamling (vinkelstykke, ventil, ventilklemme, drejestuds) / Vinkelrør (vinkelrør, ventil, ventilklemme, svivel), / Kulmakappalerakenne (kulmakappale, venttiili, venttiilin kiinnike, pyörivä liitin) / Könyök szerelvény (könyök, szelep, szelep csat, forgó csukló) / Zespół kolanka (kolanko, zawór, gniazdo zaworu, króciec do podłączenia węża) / Συγκρότημα γωνιακού εξαρτήματος (Γωνιακό εξάρτημα, Βαλβίδα, Κλιπ βαλβίδας, Στροφέας)

61294



XS* Extra Small / Erittäin pieni / Extra kicsi / Rozmiar bardzo mały / Πάρα πολύ μικρό μέγεθος

S* Small / Pieni / Kicsi / Rozmiar mały / Μικρό μέγεθος

M* Medium / Keskikokoinen / Közepes / Rozmiar średni / Μεσαίο μέγεθος

L* Large / Suuri / Nagy / Rozmiar duży / Μεγάλο μέγεθος

St* Standard / Normaali / Szabvány / Rozmiar standardowy / Κανονικό



For details, refer to the *Mirage Quattro User Guide*. Specifications may change without notice.

Mirage and Mirage Quattro are trademarks of ResMed Pty Ltd and Mirage is registered in U.S. Patent and Trademark Office. © 2012 ResMed. 618221/1 2012-01

CE
0123



618221