



## Mirage Liberty™

### FULL FACE MASK

#### FULL FACE FREEDOM

**Minimalist:** bringing minimalist design to a full face mask, and offering a sense of freedom, Mirage Liberty provides an alternative to conventional full face masks.

**Comfortable:** Mirage™ dual-wall cushion and nasal pillows with trampoline action enhance comfort around the mouth and at the nose.

**Stable:** supportive headgear keeps the mask securely in place, providing stable performance.

**Light:** one of the lightest full face masks you can find, Mirage Liberty weighs considerably less than similar masks as well as conventional full face masks.



**Dual-wall nasal pillows with trampoline action**  
adapt to different nose-to-mouth heights, angles and nasal contours

**Supportive headgear with stabilising elements**  
stays in shape when removed and provides stability when user moves during sleep

**Built-in vents**  
provide minimal noise output

**Quick-release elbow**  
allows for easy disconnection from mask

**Anti-asphyxia valve (AAV)**  
enhances safety if the airflow to the mask is impeded

**Mirage dual-wall cushion**  
provides exceptional seal and accommodates jaw movement, eliminating need for chin strap

**Set-and-forget headgear clips**  
allow convenient removal of mask

## TECHNOLOGY HIGHLIGHT

### Dual-wall nasal pillows with trampoline action

- Thin outer wall is soft and comfortable; thick inner wall ensures effective seal.
- Airflow between the two walls helps pillows adjust to nasal contours.
- Trampoline action aligns pillows to the angle of the nose and the distance between the nose and mouth.
- Pillows sit gently at the entrance to the nostrils, each moving independently of the other.



Pillows can be aligned to the nose as well as each individual nostril. The pillows can be positioned at various angles to suit individual facial characteristics.



## USER BENEFITS

### Less mask, more freedom

- Free field of vision
- Freedom to read, watch television and wear glasses
- Less contact with skin, reducing the chances of irritation
- Feels light on the face

### Comfort meets performance

- Contoured mouth cushion accommodates jaw movement and mouth leak
- Greater stability and comfort from supportive headgear with stabilising elements

### Making it easier

- Dual-wall mouth cushion and nasal pillows with trampoline action make it easy to achieve a better fit—first time, every time
- Fewer mask parts for easy assembly and disassembly

## PROFESSIONAL BENEFITS

### Fit range increased to over 90% of users

- Full face mask option for over 90% of users (based on studies of different facial types around the world)\*

### First-time fit

- First-time fit reduces callbacks and mask refitting
- Dual-wall mouth cushion and nasal pillows with trampoline action help fit users first time, every time

### Reduced inventory

- Only two product codes for the complete mask system
- Minimal mask parts

### Comprehensive user support package

- Travel and storage bag
- Mask care DVD
- Quick-fitting guide
- User guide

## PRODUCT CODES

### Mirage Liberty Full Face Mask

(three pillow sizes S, M, L packed with each mask)

#### Europe 1<sup>1</sup>

Small	61302
Large	61303

#### Europe 2<sup>2</sup>

Small	61304
Large	61305

#### Europe 3<sup>3</sup>

Small	61314
Large	61315

#### France

Small	61306
Large	61307

#### UK

Small	61308
Large	61309

#### Asia-Pacific

Small	61312
Large	61313

#### ANZ

Small	61310
Large	61311

<sup>1</sup> English, German, French, Spanish, Italian, Portuguese and Dutch

<sup>2</sup> English, Swedish, Norwegian, Danish, Finnish, Hungarian, Polish and Greek

<sup>3</sup> English, Russian, Turkish, Czech, Estonian, Slovenian, Slovakian, Bulgarian, Romanian, Latvian and Lithuanian

\* Source: internal ResMed anthropometric data