

Designed **with care**



As the leader in mask technology, ResMed has partnered with paediatric health experts to bring you Pixi™—a paediatric mask designed with children and their families in mind.

## A mask designed **specifically for children**

Uniquely constructed using the anthropometric data of children aged two and over, Pixi Paediatric Mask is not merely a scaled down adult mask.

Taking into account the bone structure, facial characteristics and skin sensitivity of children, the mask is designed to minimise pressure on soft faces.

And when combined with ResMed's latest Stellar™ ventilators or VPAP™ S or VPAP™ ST device, you have the ideal paediatric solution that delivers greater comfort and reliability for children aged two years and over weighing more than 13 kg.

Pixi Paediatric Mask is made with both the child and their carer in mind.

Every feature is carefully considered to address and improve therapy **comfort, acceptance and overall experience for everyone involved**—so carers can be free to focus on emotionally supporting their child on their journey.



**Pixi™**  
PAEDIATRIC MASK



## Technology Highlight

### NEW Headgear designed especially for children\*

**Three points of adjustment**—the mask can be adjusted at the top, side and bottom straps, providing a greater range of adjustment to cater for children of various ages.

**Out-of-sight**—the headgear is designed to sit away from the child's eyes and ears to reduce irritation and obstruction.

**Clear line of sight**—not only can the child see clearly while wearing the mask, but carers are also reassured and comforted by being able to see their child's eyes and more of their face.

*\*Designed using the specific head and facial anthropometric data of children.*

#### Soft, comfortable cushion

Single-wall, thin silicone cushion minimises pressure on children's soft faces, without compromising seal

#### Two tubing positions

Air tubing can be moved to either side of the cushion depending on your child's sleeping position

#### Spring-flex lightweight tube

stretches and bends easily, so it won't get caught or drag as your child moves around during sleep



#### Comfortable Breathaprene™ headgear materials

lightly coloured for a more reassuring, unobtrusive appearance on the child's face

#### Side-facing back strap

ergonomically designed from the fitting adult's point of view

#### Flexible stabilisers

allow the mask to conform to each child's unique facial structure and shape for the most customised fit

#### Quick-release latch

makes it easy to take the mask off, while retaining desired headgear settings

#### Emergency latch

Release the latch even faster in an emergency by pulling firmly on the lower headgear strap

## Complete System

Pixi

61031

Combine the Pixi mask with the latest Stellar or VPAP devices for a therapy system that delivers greater comfort, reliability and effective care.



Stellar



VPAP S



VPAP ST



Pixi Mask



**Pixi™**  
PAEDIATRIC MASK

Manufacturer: ResMed Ltd 1 Elizabeth Macarthur Drive Bella Vista NSW 2153 Australia. Distributed by: ResMed Corp 9001 Spectrum Center Boulevard San Diego CA 92123 USA.  
 [EC REP] ResMed (UK) Ltd 96 Milton Park Abingdon Oxfordshire OX14 4RY UK. See [www.resmed.com](http://www.resmed.com) for other ResMed locations worldwide.

Pixi is a trademark of ResMed Ltd and Velcro is a trademark of Velcro Industries B.V. © 2011 ResMed Ltd. Specifications may change without notice. 1014884/1 2011-07