

RESMED

Less mask.
More you.



The Swift™ FX for Her combines soft, lightweight design with a fit range created especially for women. It's finished in a feminine pink colour scheme, just for her.

Benefits for you ...

Designed to fit female patients

Swift FX for Her comes with XS, S and M sized pillows providing a broader fit range for female users. Once the right pillows size is selected, it requires minimal or often, no adjustment. So it saves you time, money and additional patient visits.

Interchangeable with Swift FX™ Bella

Swift FX for Her is fully interchangeable with Swift FX Bella (both use the same pillows and tubing) for an increased level of patient compliance.

Easier therapy acceptance

Swift FX for Her is not intimidating, helping users accept and continue therapy. With very few parts, Swift FX for Her is easy to clean and assemble too.

Gentle therapy, stable seal

Flexible pillows with trampoline action allow movement in multiple directions while maintaining a seal.

... and your patients

Great hair management

The length of the backstrap is shorter than the standard Swift FX. It's designed so that women can choose to wear their hair over or under the strap to accommodate a variety of hairstyles.

Interchangeable with Swift FX Bella

It is fully interchangeable with Swift FX Bella (both have the same pillows and tubing) which provides users with another headgear option.

Feminine colour scheme

Low-profile backstrap with Velcro™ fasteners and softwraps in feminine pink.

Quiet performance

Innovative vent design diffuses airflow for quiet mask performance—only 25 dBA.*

* Declared dual numbered noise emission values (with uncertainty 3dBA) in accordance with ISO 4871.e 1



Swift™ FX for Her
NASAL PILLOWS SYSTEM



There's nothing hard about it...

The Swift FX for Her is so soft, light and simple, it feels like it's barely there.

Minimal contact on the face and unobtrusive design make therapy less intimidating and easy to accept.

Pillows for her (XS, S and M) come packaged with the mask system to fit most female patients.

Dual-wall nasal pillows seal on contact for a gentle, comfortable seal, offer a wide fit range and disperse air gently to the nose

Soft, narrow cushion with integrated flexible chamber makes side sleeping comfortable and helps maintain seal during movement

Ultra-smooth 360° rotating elbow rotates with easy, fluid movement for flexibility in tube positions

'Lock and key' for easy attaching and detaching of headgear

Product Codes

Europe 1	61545
Europe 2	61546
Europe 3	61547

Europe 1: English (UK), French, German and Italian (Switzerland)

Europe 2: English, Danish, Finnish, Norwegian and Swedish

Europe 3: English, Czech, Dutch (Belgium), French (Belgium), German (Belgium), Greek, Hungarian, Italian, Polish, Portuguese, Russian, Slovak, Spanish and Turkish

Spring-flex lightweight tube stretches and bends uniquely for increased stability without the use of a tube retainer

Quiet vents diffuse and direct airflow quietly away from the user and bed partner

Simple top buckle has clearly marked notches for easy adjustment

Soft silicone headgear embraces the face gently and requires minimal adjustment

Low-profile backstrap with Velcro fasteners is available in a small size, in a feminine pink colour scheme

Soft wraps enhance comfort at the cheeks. Soft wraps available in pink to colour coordinate with the backstrap

Size and name indicators make it easy to identify the size and the correct orientation of the mask

Flexible pillows with trampoline action accommodate movement in multiple directions while maintaining seal

Swift™ FX for Her
NASAL PILLOWS SYSTEM



ResMed Ltd
1 Elizabeth Macarthur Drive
Bella Vista NSW 2153 Australia

Distributed by: ResMed Corp 9001 Spectrum Center Boulevard San Diego CA 92123 USA. **EC REP:** ResMed (UK) Ltd 96 Milton Park Abingdon Oxfordshire OX14 4RY UK. See www.resmed.com for other ResMed locations worldwide. Swift is a trademark of ResMed Ltd and is registered in U.S. Patent and Trademark Office. Velcro is a trademark of Velcro Industries B.V. © 2012 ResMed Ltd. Specifications may change without notice. 1015224/2 2012-02