

Ultra Mirage™

NON-VENTED FULL FACE MASK

The stable choice for noninvasive ventilation


ResMed's Ultra Mirage Non-vented (NV) Full Face Mask provides superior seal, stability and comfort for long-term noninvasive ventilation in the home as well as the hospital environment.

The Ultra Mirage NV Full Face Mask provides all the advantages of the original Mirage™ Non-vented Full Face Mask.

It also offers new features and design improvements to enhance therapy, ease-of-use and comfort.

The unique adjustable forehead support enables fitting to suit each individual patient.

The most effective position can be selected while the patient is wearing the mask. This helps maintain seal and improves overall stability.



TECHNOLOGY HIGHLIGHT

Adjustable forehead support—provides an effective seal, even when the patient moves or opens their mouth

- The forehead support, with four adjustable positions, provides an ideal seal and enhances comfort.
- Selecting the most effective position helps maintain seal even during minor patient movements, such as the patient opening their mouth.
- The flexible forehead pads conform to the contour of the patient's forehead. The pads distribute pressure evenly over a larger surface area, reducing pressure on the patient's face.



Patient benefits

Proven technology, trusted comfort

- Mirage dual-wall cushion technology reduces pressure on the nose, improving seal and comfort.
- 360° rotating elbow provides extra stability and convenience.

Reassurance that prescribed therapy is being delivered

- Superior seal and stability offer peace of mind for patients.

User convenience

- Quick-release headgear clips enable quick and easy mask fitting and removal.
- Once set up, the headgear doesn't have to be readjusted.

Professional benefits

Confidence in effective therapy delivery

- As the mask covers both the nose and mouth, it provides effective therapy even in the presence of mouth breathing.

Tried and trusted with proven quality and technology

- ResMed's patented Mirage dual-wall cushion technology encourages compliance by enhancing seal and comfort.

Quick-release headgear clips

- Allow easy removal of the mask from the headgear without disturbing the patient.
- Four clips, instead of two, make fitting easier and quicker.

Efficient mask size selection

- 3D fitting templates help achieve first-time fit.
- Use of templates can also reduce mask fitting time.

PRODUCT CODES

Ultra Mirage Non-vented Full Face Mask

Europe (Group 1¹), ANZ, Asia Pacific

Small Shallow	60642
Small Standard	60643
Medium Shallow	60644
Medium Standard	60645
Large Shallow	60646
Large Standard	60647

Europe (Group 2²)

Small Standard	60648
Small Shallow	60649
Medium Standard	60650
Medium Shallow	60651
Large Standard	60652
Large Shallow	60653

European user information available in the following languages:
¹ English, German, French, Spanish, Italian, Portuguese, and Dutch
² English, Swedish, Norwegian, Danish, Finnish