

# Zoznam kompatibilných másk a prístrojov

Slovensky

Tento dokument špecifikuje masky, ktoré sú kompatibilné s prístrojom pozitívneho pretlaku ResMed a prípadné potrebné nastavenia masky. Pre podrobné informácie o správnom používaní týchto prístrojov pozrite špecifický návod na používanie.

Poznámka: Nie všetky masky sú k dispozícii vo všetkých oblastiach.

## Celotvárové masky

Prístroj	AcuCare F1-4	AirFit™ F30i <sup>1</sup>	AirFit™ F30	AirTouch™ F20 / AirTouch™ F20 for Her / AirFit™ F20 / AirFit™ F20 for Her	AirFit™ F10 / AirFit™ F10 for Her	Mirage™	Mirage Series 2	Ultra Mirage™	Mirage Quattro™	Quattro™ FX / Quattro™ FX for Her / Quattro™ Air / Quattro™ Air for Her	Mirage Liberty™	Hospital Full Face
AirMini™	X	X	✓	✓	X	X	X	X	X	X	X	X
AirStart™ 10 / AirSense™ 10	N/S	Full Face	Full Face	Full Face	Full Face	N/S	N/S	Full Face	Full Face	Full Face	Full Face	Full Face
AirCurve™ 10 VAuto / AirCurve™ 10 S / AirCurve™ 10 ST	N/S	Full Face	Full Face	Full Face	Full Face	N/S	N/S	Full Face	Full Face	Full Face	Full Face	Full Face
AirCurve™ 10 CS Pacewave	N/S	Full Face	Full Face	Full Face	Full Face	N/S	N/S	Full Face	Full Face	Full Face	Full Face	Full Face
AirCurve™ 10 CS-A Pacewave	N/S	N/S	Full Face	Full Face	Full Face	N/S	N/S	Full Face	Full Face	Full Face	Full Face	Full Face
Lumis™	N/S	N/S	Full Face <sup>3</sup>	Full Face	Full Face	N/S	N/S	Full Face	Full Face	Full Face	Full Face	Full Face
S9 AutoSet™ / S9 Elite™	N/S	Full Face	N/S	Full Face	Full Face	N/S	N/S	Full Face	Full Face	Full Face	Full Face	Full Face
S9 Escape™	N/S	Full Face	N/S	Full Face	Full Face	N/S	N/S	Full Face	Full Face	Full Face	Full Face	Full Face
AutoSet CS™ 2	N/S	N/S	N/S	N/S	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE <sup>1</sup>	FULL FACE
AutoSet CS™ (S9)	N/S	Full Face	N/S	Full Face	Full Face	N/S	N/S	Full Face	Full Face	Full Face	Full Face	Full Face
AutoSet CS™ (S9 PaceWave)	N/S	Full Face	N/S	Full Face	Full Face	N/S	N/S	Full Face	Full Face	Full Face	Full Face	Full Face
AutoSet CS-A™ (S9 PaceWave)	N/S	N/S	N/S	Full Face	Full Face	N/S	N/S	Full Face	Full Face	Full Face	Full Face	Full Face
C-Series Tango™	N/S	N/S	N/S	N/S	Nie sú potrebné žiadne nastavenia. Poradte sa s Vaším lekárom pred prechodom k inému typu masky.							
S7™ Lightweight	N/S	N/S	N/S	N/S	Nie sú potrebné žiadne nastavenia. Poradte sa s Vaším lekárom pred prechodom k inému typu masky.							
S7 Elite™ / S7 AutoSet Spirit	N/S	N/S	N/S	N/S	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
S7 AutoSet Respond™	N/S	N/S	N/S	N/S	MIRAGE	STANDARD	ULTRA	MIRAGE	MIRAGE	MIRAGE	MIR FULL	MIR FULL
S8™ Auto 25	N/S	N/S	N/S	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
S8 AutoSet Spirit™	N/S	N/S	N/S	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
S8 AutoSet Vantage™	N/S	N/S	N/S	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
S8 Escape™	N/S	N/S	N/S	N/S	Vyžaduje si tlakovú kalibráciu pri 12 cm H2O. Poradte sa so svojím lekárom a pozrite si používateľskú príručku.							
S8 Escape™ II	N/S	N/S	N/S	Vyžaduje si tlakovú kalibráciu pri 12 cm H2O. Poradte sa so svojím lekárom a pozrite si používateľskú príručku.								
S8 Elite™ / S8 Elite™ II	N/S	N/S	N/S	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
S8 AutoSet Spirit™ II	N/S	N/S	N/S	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL

Prístroj	AcuCare F1-4	AirFit™ F30i <sup>1</sup>	AirFit™ F30	AirTouch™ F20 / AirTouch™ F20 for Her AirFit™ F20 / AirFit™ F20 for Her	AirFit™ F10 / AirFit™ F10 for Her	Mirage™	Mirage Series 2	Ultra Mirage™	Mirage Quattro™	Quattro™ FX / Quattro™ FX for Her Quattro™ Air / Quattro™ Air for Her	Mirage Liberty™	Hospital Full Face
S8 Compact™	N/S	N/S	N/S	N/S		Nie sú potrebné žiadne nastavenia. Poradte sa s Vaším lekárom pred prechodom k inému typu masky.						
VPAP™ III / VPAP™ III ST	N/S	N/S	N/S	N/S	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
VPAP™ III ST-A	N/S	N/S	N/S	N/S	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
VPAP™ IV / VPAP™ IV ST	N/S	N/S	N/S	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
VPAP™ Adapt SV	N/S	N/S	N/S	N/S	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE <sup>1</sup>	FULL FACE
VPAP™ Adapt (S9)	N/S	Full Face	N/S	Full Face	Full Face	N/S	N/S	Full Face	Full Face	Full Face	Full Face	Full Face
VPAP™ Adapt (S9 PaceWave)	N/S	Full Face	N/S	Full Face	Full Face	N/S	N/S	Full Face	Full Face	Full Face	Full Face	Full Face
VPAP™ Malibu	N/S	N/S	N/S	N/S	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE
VPAP™ Auto 25	N/S	N/S	N/S	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
Auto 25 (S9)	N/S	Full Face	N/S	Full Face	Full Face	N/S	N/S	Full Face	Full Face	Full Face	Full Face	Full Face
VPAP™ S / VPAP™ ST	N/S	N/S	N/S	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
VPAP™ S (S9) / VPAP™ ST (S9)	N/S	Full Face	N/S	Full Face	Full Face	N/S	N/S	Full Face	Full Face	Full Face	Full Face	Full Face
VPAP™ ST-A (S9)	N/S	N/S	N/S	Full Face	Full Face	N/S	N/S	Full Face	Full Face	Full Face	Full Face	Full Face
VPAP™ ST (S9 iVAPS) / VPAP™ ST-A (S9 iVAPS)	N/S	N/S	N/S	Full Face	Full Face	N/S	N/S	Full Face	Full Face	Full Face	Full Face	Full Face
S9 VPAP™ Tx	N/S	Full Face <sup>2</sup>	N/S	Full Face	Full Face	N/S	N/S	Full Face	Full Face	Full Face	Full Face	Full Face
VPAP™ Tx	N/S	N/S	N/S	N/S	FULL FACE	N/S	N/S	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE
Stellar™	N/S	N/S	N/S	FULL FACE	FULL FACE	N/S	N/S	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE
VS Ultra™ / VS Integra™	✓	N/S	N/S	N/S	✓	N/S	N/S	✓	✓	N/S	✓	✓
VS III™	N/S	N/S	N/S	N/S	✓	N/S	N/S	✓	✓	N/S	✓	✓
Astral™	FULL FACE	N/S	N/S	FULL FACE	N/S	N/S	N/S	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE

<sup>1</sup> Če to masko uporabljate z nekaterimi napravami CPAP ali napravami za dvostopenjsko zdravljenje, SmartStop morda ne bo deloval, kot je opisano.

<sup>2</sup> Brez načinov PAC in iVAPS. Združljiva samo s tlakom v razponu 4–25 cm H2O.

<sup>3</sup> Účinnosť liečby môže byť ovplyvnená v oblastiach, ako je presnosť kontroly AutoEPAP, a to najmä pri vyšších tlakoch.

Upozornenie: Je potrebné overiť účinnosť ventilácie a alarmov vrátane po každej zmene nastavenia ventilácie a alarmov, konfigurácii obvodu a prídruženej liečbe (napríklad nebulizácia, prívod kyslíka)

N/S = maska nie je schválená na používanie s týmto prístrojom.

✓ = Kompatibilná X = Nekompatibilná

## Nosné masky 1

Prístroj	AirFit™ N30 for AirMini	AirFit™ N30i <sup>1</sup> / AirFit™ N30	AirFit™ N20 Classic	AirFit™ N20 / AirFit™ N20 for Her	AirTouch™ N20 / AirTouch™ N20 for Her	AirFit™ N10 / AirFit™ N10 for Her	Swift™ FX Nano / Swift™ FX Nano for Her	Mirage™ FX / Mirage FX™ for Her	Mirage Activa™ LT	Mirage™ SoftGel	Hospital Nasal
AirMini™	✓	X	✓	✓	✓	X	X	X	X	X	X
AirStart™ 10 / AirSense™ 10	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
AirCurve™ 10 VAuto / AirCurve™ 10 S / AirCurve™ 10 ST	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
AirCurve™ 10 CS Pacewave	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
AirCurve™ 10 CS-A Pacewave	X	N/S	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
Lumis™	X	N/S	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
S9 AutoSet™ / S9 Elite™	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal

Prístroj	AirFit™ N30 for AirMini	AirFit™ N30i <sup>1</sup> / AirFit™ N30	AirFit™ N20 Classic	AirFit™ N20 / AirFit™ N20 for Her	AirTouch™ N20 / AirTouch™ N20 for Her	AirFit™ N10 / AirFit™ N10 for Her	Swift™ FX Nano / Swift™ FX Nano for Her	Mirage™ FX / Mirage FX™ for Her	Mirage Activa™ LT	Mirage™ SoftGel	Hospital Nasal
S9 Escape™	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
AutoSet CS™ 2	X	N/S	N/S	N/S	N/S	N/S	N/S	VISTA <sup>1</sup>	ACTIVA	ACTIVA	VISTA
AutoSet CS™ (S9)	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
AutoSet CS™ (S9 PaceWave)	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
AutoSet CS-A™ (S9 PaceWave)	X	N/S	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
C-Series Tango™	X	N/S	N/S	N/S	N/S	N/S	N/S	Niesú potrebné žiadne nastavenia. Poradte sa s Vaším lekárom pred prechodom kinémutypumasky.			
S7™ Lightweight	X	N/S	N/S	N/S	N/S	N/S	N/S	Niesú potrebné žiadne nastavenia. Poradte sa s Vaším lekárom pred prechodom kinémutypumasky.			
S7 Elite™ / S7 AutoSet Spirit	X	N/S	N/S	N/S	N/S	N/S	N/S	STANDARD	STANDARD	STANDARD	STANDARD
S7 AutoSet Respond™	X	N/S	N/S	N/S	N/S	N/S	N/S	STANDARD	STANDARD	STANDARD	STANDARD
S8™ Auto 25	X	N/S	N/S	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	ACTIVA	ACTIVA	STANDARD
S8 AutoSet Spirit™	X	N/S	N/S	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	ACTIVA	ACTIVA	STANDARD
S8 AutoSet Vantage™	X	N/S	N/S	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	ACTIVA	ACTIVA	STANDARD
S8 Escape™	X	N/S	N/S	N/S	N/S	Vyžaduje si tlakovú kalibráciu pri 12 cm H2O. Poradte sa so svojím lekárom a pozrite si používateľskú príručku.					
S8 Escape™ II	X	N/S	N/S	Vyžaduje si tlakovú kalibráciu pri 12 cm H2O. Poradte sa so svojím lekárom a pozrite si používateľskú príručku.							
S8 Elite™ / S8 Elite™ II	X	N/S	N/S	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	ACTIVA	ACTIVA	STANDARD
S8 AutoSet Spirit™ II	X	N/S	N/S	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	ACTIVA	ACTIVA	STANDARD
S8 Compact™	X	N/S	N/S	N/S	N/S	Nie sú potrebné žiadne nastavenia. Poradte sa s Vaším lekárom pred prechodom k inému typu masky.					
VPAP™ III / VPAP™ III ST	X	N/S	N/S	N/S	N/S	MIRAGE	MIRAGE	STANDARD	STANDARD	STANDARD	N/S
VPAP™ III ST-A	X	N/S	N/S	N/S	N/S	MIRAGE	MIRAGE	STANDARD	STANDARD	STANDARD	N/S
VPAP™ IV / VPAP™ IV ST	X	N/S	N/S	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	ACTIVA	ACTIVA	STANDARD
VPAP™ Adapt SV	X	N/S	N/S	N/S	N/S	N/S	N/S	VISTA <sup>1</sup>	ACTIVA	ACTIVA	VISTA
VPAP™ Adapt (S9)	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
VPAP™ Adapt (S9 PaceWave)	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
VPAP™ Malibu	X	N/S	N/S	N/S	N/S	SWIFT	SWIFT	VISTA <sup>1</sup>	N/S	N/S	MERIDIAN
VPAP™ Auto 25	X	N/S	N/S	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	ACTIVA	ACTIVA	STANDARD
Auto 25 (S9)	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
VPAP™ S / VPAP™ ST	X	N/S	N/S	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	ACTIVA	ACTIVA	STANDARD
VPAP™ S (S9) / VPAP™ ST (S9)	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
VPAP™ ST-A (S9)	X	N/S	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
VPAP™ ST (S9 iVAPS) / VPAP™ ST-A (S9 iVAPS)	X	N/S	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
S9 VPAP™ Tx	X	Pillows <sup>2</sup>	N/S	Pillows	N/S	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
VPAP™ Tx	X	N/S	N/S	N/S	N/S	PILLOWS	PILLOWS	NASAL	NASAL	NASAL	NASAL
Stellar™	X	N/S	PILLOWS	PILLOWS	PILLOWS	PILLOWS <sup>1</sup>	PILLOWS <sup>1</sup>	NASAL	NASAL	NASAL	NASAL
VS Ultra™ / VS Integra™	X	N/S	N/S	N/S	N/S	N/S	N/S	N/S	N/S	N/S	N/S
VS III™	X	N/S	N/S	N/S	N/S	N/S	N/S	N/S	✓	N/S	✓
Astral™	X	N/S	Pillows	Pillows	Pillows	N/S	N/S	Nasal	N/S	Nasal	Nasal

<sup>1</sup> Če to masko uporabljate z nekaterimi napravami CPAP ali napravami za dvostopenjsko zdravljenje, SmartStop morda ne bo deloval, kot je opisano.

<sup>2</sup> Brez načinov PAC in iVAPS. Zdržljiva samo s tlakom v razponu 4–25 cm H2O.

N/S = maska nie je schválená na používanie s týmto prístrojom.

## Nosné masky 2

Prístroj	Pixi™	Mirage Kidsta™	Mirage Micro™	Mirage Vista™	Ultra Mirage II	Mirage	Mirage Activa™	Modular	Papillon™	Meridian
AirMini™	X	X	X	X	X	X	X	X	X	X
AirStart™ 10 / AirSense™ 10	N/S	Nasal	Nasal	Nasal	Nasal	N/S	Nasal	N/S	N/S	Nasal
AirCurve™ 10 VAuto / AirCurve™ 10 S / AirCurve™ 10 ST	N/S	Nasal	Nasal	Nasal	Nasal	N/S	Nasal	N/S	N/S	Nasal
AirCurve™ 10 CS Pacewave	N/S	Nasal	Nasal	Nasal	Nasal	N/S	Nasal	N/S	N/S	Nasal
AirCurve™ 10 CS-A Pacewave	N/S	Nasal	Nasal	Nasal	Nasal	N/S	Nasal	N/S	N/S	Nasal
Lumis™	N/S	Nasal	Nasal	Nasal	Nasal	N/S	Nasal	N/S	N/S	Nasal
S9 AutoSet™ / S9 Elite™	N/S	N/S	Nasal	Nasal	Nasal Ultra	N/S	Nasal	Nasal	Nasal	Nasal
S9 Escape™	N/S	N/S	Nasal	Nasal	Nasal Ultra	N/S	Nasal	Nasal	Nasal	Nasal
AutoSet CS™ 2	N/S	N/S	VISTA¹	VISTA	ULTRA	N/S	ACTIVA	N/S	PAPILLON	VISTA
AutoSet CS™ (S9)	N/S	N/S	Nasal	Nasal	Nasal Ultra	N/S	Nasal	Nasal	Nasal	Nasal
AutoSet CS™ (S9 PaceWave) / AutoSet CS-A™ (S9 PaceWave)	N/S	N/S	Nasal	Nasal	Nasal Ultra	N/S	Nasal	Nasal	Nasal	Nasal
C-Series Tango™	N/S	N/S	N/S	Nie sú potrebné žiadne nastavenia. Poradte sa s Vaším lekárom pred prechodom k inému typu masky.						N/S
S7™ Lightweight	N/S	N/S	N/S	Nie sú potrebné žiadne nastavenia. Poradte sa s Vaším lekárom pred prechodom k inému typu masky.						N/S
S7 Elite™ / S7 AutoSet Spirit	N/S	N/S	STANDARD	STANDARD	ULTRA	MIRAGE	STANDARD	STANDARD	STANDARD	N/S
S7 AutoSet Respond™	N/S	N/S	STANDARD	STANDARD	ULTRA	MIRAGE	STANDARD	STANDARD	STANDARD	N/S
S8™ Auto 25	N/S	N/S	STANDARD	STANDARD	ULTRA	MIRAGE	ACTIVA	N/S	N/S	STANDARD
S8 AutoSet Spirit™	N/S	N/S	STANDARD	STANDARD	ULTRA	MIRAGE	ACTIVA	STANDARD	STANDARD	STANDARD
S8 AutoSet Vantage™	N/S	N/S	STANDARD	STANDARD	ULTRA	MIRAGE	ACTIVA	STANDARD	STANDARD	STANDARD
S8 Escape™	N/S	N/S	Vyžaduje si tlakovú kalibráciu pri 12 cm H2O. Poradte sa so svojím lekárom a pozrite si používateľskú príručku.							
S8 Escape™ II	N/S	N/S	Vyžaduje si tlakovú kalibráciu pri 12 cm H2O. Poradte sa so svojím lekárom a pozrite si používateľskú príručku.							
S8 Elite™ / S8 Elite™ II	N/S	N/S	STANDARD	STANDARD	ULTRA	MIRAGE	ACTIVA	STANDARD	STANDARD	STANDARD
S8 AutoSet Spirit™ II	N/S	N/S	STANDARD	STANDARD	ULTRA	MIRAGE	ACTIVA	STANDARD	STANDARD	STANDARD
S8 Compact™	N/S	N/S	Nie sú potrebné žiadne nastavenia. Poradte sa s Vaším lekárom pred prechodom k inému typu masky.							
VPAP™ III / VPAP™ III ST	N/S	N/S	STANDARD	STANDARD	ULTRA	MIRAGE	STANDARD	STANDARD	N/S	N/S
VPAP™ III ST-A	N/S	STANDARD	STANDARD	STANDARD	ULTRA	MIRAGE	STANDARD	STANDARD	N/S	N/S
VPAP™ IV / VPAP™ IV ST	N/S	N/S	STANDARD	STANDARD	ULTRA	MIRAGE	ACTIVA	N/S	N/S	STANDARD
VPAP™ Adapt SV	N/S	N/S	VISTA¹	VISTA	ULTRA	N/S	ACTIVA	N/S	N/S	VISTA
VPAP™ Adapt (S9)	N/S	N/S	Nasal	Nasal	Nasal Ultra	N/S	Nasal	Nasal	Nasal	Nasal
VPAP™ Adapt (S9 PaceWave)	N/S	N/S	Nasal	Nasal	Nasal Ultra	N/S	Nasal	Nasal	Nasal	Nasal
VPAP™ Malibu	N/S	N/S	MERIDIAN	VISTA	ULTRA	N/S	N/S	N/S	N/S	MERIDIAN
VPAP™ Auto 25	N/S	N/S	STANDARD	STANDARD	ULTRA	MIRAGE	ACTIVA	N/S	N/S	STANDARD
Auto 25 (S9)	N/S	N/S	Nasal	Nasal	Nasal Ultra	N/S	Nasal	Nasal	Nasal	STANDARD
VPAP™ S / VPAP™ ST	N/S	N/S	STANDARD	STANDARD	ULTRA	MIRAGE	ACTIVA	N/S	N/S	STANDARD
VPAP™ S (S9) / VPAP™ ST (S9) / VPAP™ ST-A (S9)	Pediatric	N/S	Nasal	Nasal	Nasal Ultra	N/S	Nasal	Nasal	Nasal	Nasal

Prístroj	Pixi™	Mirage Kidsta™	Mirage Micro™	Mirage Vista™	Ultra Mirage II	Mirage	Mirage Activa™	Modular	Papillon™	Meridian
VPAP™ ST (S9 iVAPS) / VPAP™ ST-A (S9 iVAPS)	Pediatric	N/S	Nasal	Nasal	Nasal Ultra	N/S	Nasal	Nasal	Nasal	Nasal
S9 VPAP™ Tx	Pediatric	N/S	Nasal	Nasal	Nasal Ultra	N/S	Nasal	N/S	N/S	Nasal
VPAP™ Tx	N/S	NASAL	NASAL	NASAL	NASAL-ULTRA	N/S	NASAL	N/S	N/S	NASAL
Stellar™	Pediatric	NASAL	NASAL	NASAL	ULTRA	N/S	NASAL	N/S	N/S	N/S
VS Ultra™ / VS Integra™	N/S	N/S	N/S	N/S	N/S	N/S	N/S	N/S	N/S	N/S
VS III™	N/S	✓	✓	✓	✓	N/S	✓	N/S	N/S	N/S
Astral™	Pediatric	Nasal	Nasal	Nasal	Nasal Ultra	N/S	N/S	N/S	N/S	N/S

N/S = maska nie je schválená na používanie s týmto prístrojom.

✓ = Kompatibilná X = Nekompatibilná

## Nosné vankúšiky

Prístroj	AirFit™ P30i <sup>1</sup>	AirFit™ P10 For AirMin™	AirFit™ P10 / AirFit™ P10 for Her	Mirage Swift™	Mirage Swift™ II	Swift™ FX / Swift™ FX for Her / Swift™ FX Bella / Swift™ FX Bella Gray	Swift™ LT	Swift™ LT for Her / Swift™ LT-F		
AirMini™	X	✓	X	X	X	X	X	X		
AirStart™ 10 / AirSense™ 10	Pillows	X	Pillows <sup>2</sup>	N/S	Pillows	Pillows	Pillows	Pillows		
AirCurve™ 10 VAuto / AirCurve™ 10 S / AirCurve™ 10 ST	Pillows	X	Pillows <sup>2</sup>	N/S	Pillows	Pillows	Pillows	Pillows		
AirCurve™ 10 CS Pacewave	Pillows	X	Pillows <sup>2</sup>	N/S	Pillows	Pillows	Pillows	Pillows		
AirCurve™ 10 CS-A Pacewave	N/S	X	Pillows <sup>2</sup>	N/S	Pillows	Pillows	Pillows	Pillows		
Lumis™	N/S	X	Pillows <sup>2</sup>	N/S	Pillows	Pillows	Pillows	Pillows		
S9 AutoSet™ / S9 Elite™	Pillows	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows		
S9 Escape™	Pillows	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows		
AutoSet CS™ 2	N/S	X	N/S	N/S	N/S	N/S	N/S	N/S		
AutoSet CS™ (S9)	Pillows	X	Pillows <sup>2</sup>	Pillows	Pillows	Pillows	Pillows	Pillows		
AutoSet CS™ (S9 PaceWave)	Pillows	X	Pillows <sup>2</sup>	Pillows	Pillows	Pillows	Pillows	Pillows		
AutoSet CS-A™ (S9 PaceWave)	N/S	X	Pillows <sup>2</sup>	Pillows	Pillows	Pillows	Pillows	Pillows		
C-Series Tango™	N/S	X		Nie sú potrebné žiadne nastavenia. Poradte sa s Vaším lekárom pred prechodom k inému typu masky.						
S7™ Lightweight	N/S	X		Nie sú potrebné žiadne nastavenia. Poradte sa s Vaším lekárom pred prechodom k inému typu masky.						
S7 Elite™ / S7 AutoSet Spirit	N/S	X	MIRAGE	MIRAGE	MIRAGE	MIRAGE	MIRAGE	MIRAGE		
S7 AutoSet Respond™	N/S	X	MIRAGE	MIRAGE	MIRAGE	MIRAGE	MIRAGE	MIRAGE		
S8™ Auto 25	N/S	X	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT		
S8 AutoSet Spirit™	N/S	X	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT		
S8 AutoSet Vantage™	N/S	X	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT		
S8 Escape™	N/S	X		Vyžaduje si tlakovú kalibráciu pri 12 cm H2O. Poradte sa so svojím lekárom a pozrite si používateľskú príručku.						
S8 Escape™ II	N/S	X		Vyžaduje si tlakovú kalibráciu pri 12 cm H2O. Poradte sa so svojím lekárom a pozrite si používateľskú príručku.						
S8 Elite™ / S8 Elite™ II	N/S	X	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT		
S8 AutoSet Spirit™ II	N/S	X	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT		
S8 Compact™	N/S	X		Nie sú potrebné žiadne nastavenia. Poradte sa s Vaším lekárom pred prechodom k inému typu masky.						
VPAP™ III / VPAP™ III ST	N/S	X	MIRAGE <sup>2</sup>	MIRAGE	MIRAGE	MIRAGE	MIRAGE	MIRAGE		

Prístroj	AirFit™ P30i <sup>1</sup>	AirFit™ P10 For AirMini™	AirFit™ P10 / AirFit™ P10 for Her	Mirage Swift™	Mirage Swift™ II	Swift™ FX / Swift™ FX for Her / Swift™ FX Bella / Swift™ FX Bella Gray	Swift™ LT	Swift™ LT for Her / Swift™ LT-F
VPAP™ III ST-A	N/S	X	MIRAGE <sup>2</sup>	MIRAGE	MIRAGE	MIRAGE	MIRAGE	MIRAGE
VPAP™ IV / VPAP™ IV ST	N/S	X	SWIFT <sup>2</sup>	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT
VPAP™ Adapt SV	N/S	X	N/S	N/S	N/S	N/S	N/S	N/S
VPAP™ Adapt (S9)	Pillows	X	Pillows <sup>2</sup>	Pillows	Pillows	Pillows	Pillows	Pillows
VPAP™ Adapt (S9 PaceWave)	Pillows	X	Pillows <sup>2</sup>	Pillows	Pillows	Pillows	Pillows	Pillows
VPAP™ Malibu	N/S	X	SWIFT <sup>2</sup>	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT
VPAP™ Auto 25	N/S	X	SWIFT <sup>2</sup>	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT
Auto 25 (S9)	Pillows	X	Pillows <sup>2</sup>	Pillows	Pillows	Pillows	Pillows	Pillows
VPAP™ S / VPAP™ ST	N/S	X	SWIFT <sup>2</sup>	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT
VPAP™ S (S9) / VPAP™ ST (S9)	Pillows	X	Pillows <sup>2</sup>	Pillows	Pillows	Pillows	Pillows	Pillows
VPAP™ ST-A (S9)	N/S	X	Pillows <sup>2</sup>	Pillows	Pillows	Pillows	Pillows	Pillows
VPAP™ ST (S9 iVAPS) / VPAP™ ST-A (S9 iVAPS)	N/S	X	Pillows <sup>2</sup>	Pillows	Pillows	Pillows	Pillows	Pillows
S9 VPAP™ Tx	Pillows <sup>3</sup>	X	Pillows <sup>2</sup>	N/S	Pillows	Pillows	Pillows	Pillows
VPAP™ Tx	N/S	X	PILLOWS <sup>2</sup>	PILLOWS	PILLOWS	PILLOWS	PILLOWS	PILLOWS
Stellar™	N/S	X	PILLOWS <sup>1, 2</sup>	PILLOWS <sup>1</sup>	PILLOWS <sup>1</sup>	PILLOWS <sup>1</sup>	PILLOWS <sup>1</sup>	PILLOWS <sup>1</sup>
VS Ultra™ / VS Integra™	N/S	X	N/S	N/S	N/S	N/S	N/S	N/S
VS III™	N/S	X	N/S	N/S	N/S	N/S	N/S	N/S
Astral™	N/S	X	N/S	Pillows	Pillows	Pillows	Pillows	Pillows

<sup>1</sup> Če to masko uporabljate z nekaterimi napravami CPAP ali napravami za dvostopenjsko zdravljenje, SmartStop morda ne bo deloval, kot je opisano.

<sup>2</sup> Presnosta respiracijskega volumena / minutovolumena dihanja se pri uporabi maske AirFit P10 v bivelel režimih lahko zniža, a to najmanj pri višji tlakovni podpori.

<sup>3</sup> Brez načinov PAC in iVAPS. Združljiva samo s tlakom v razponu 4–25 cm H<sub>2</sub>O.

N/S = maska ni je odobrena za uporabo s tem priročnikom.

✓ = Kompatibilna X = Nekompatibilna