



ResMed

LifeChoice **Activox**[™] 4L

Portable oxygen concentrator

Effective, portable oxygen therapy

The LifeChoice Activox 4L portable oxygen concentrator (POC) is a lightweight, portable, user-friendly device that delivers the right dose of oxygen exactly when it's needed.

Packed with responsive technologies, the LifeChoice Activox 4L is firmly focused on meeting your patients' needs and offering them the freedom, comfort and confidence that comes with portable therapeutic support.



Comfortable breathing on the go

The LifeChoice Activox 4L device is equipped with responsive algorithms, variable settings and relevant features that meet your patients' therapeutic requirements. Adjustable, practical carry bags and multiple charging options give them greater freedom and flexibility to live life as they choose.

Key features

Light and portable

- Lightweight at 2.2 kg (4.8 lbs)
- User-friendly carry case with straps converts easily between four carrying options

Powered to keep up

- Up to 10 hours of use with the powerful internal battery*
- An extra 4 hours of use with the optional external battery*
- AC/DC adapter delivers quick, efficient charging and operation both at home and on the go

Simple, easy-to-use interface

- Easy-to-read transfective LCD display screen
- Control panel with day-to-day information, error messages and alarms

Efficient oxygen delivery

- PULSE-WAVE™ technology effectively delivers the required dose of oxygen while minimising waste
- Auto Mode senses changes in your patient's activity level and adapts to their breathing requirements
- Auto Dose, a feature of Auto Mode, provides a 'safety' dose of oxygen

Technical specifications

Oxygen delivery	PULSE-WAVE Settings 1-4
Oxygen concentration	90% ± 3%
Weight	2.2 kg (4.8 lbs)
Size (W x H x D)	230 mm x 200 mm x 112 mm
FAA approved	•

Battery runtimes

Setting	1	2	3	4
Internal battery runtime**	10hr 15min	8hr 15min	5hr	4hr
External battery runtime**	4hr 45min	3hr 45min	2hr 45min	2hr 25min
Total runtime**	15hr	12hr	7hr 45min	6hr 25min
Charge time**	4 hours internal - 2 hours external			

Product codes

A Activox P4L (package)	AOX-P4L-EU
B Carry case	123-102
C Shoulder strap kit	XYC414
D DC charger	XYC104
E AC charger EU power cord	XYC089 XYC494
F Accessory bag	XYC340
G Adult nasal cannula (sold separately)	XYC262
H External battery (sold separately)	200122



* In setting 1

** Hours are approximate. Runtimes are based on breaths per minute.