

Resmed CPAP supplies & replacement parts

☑ Replacement tips for mask parts. What to look for:

- You have to tighten the straps often
- The mask has lost its spring
- Your therapy feels less effective than usual
- The silicone loses its grip and remains greasy, even after cleaning
- Mask seal is deteriorating or becoming more challenging

📄 Insurance coverage tip

Before ordering any supplies or replacements, we recommend contacting your insurance provider to make sure the parts you need are covered under your plan.



Replacement tips for mask parts

These replacement schedule guidelines are based on the current U.S. Medicare replacement schedule.*

☑ **Twice per month:** Nasal mask cushions, Nasal pillows mask cushions, Oral nasal mask cushions

☑ **Every 3 months:** Mask frame systems

Frame systems are in close proximity to your face and should be replaced regularly to keep therapy effective.

☑ **Every month:** Full face mask cushions

☑ **Every 6 months:** Mask headgear

Headgear can become stretched out and lose elasticity with regular use, which can lead to overtightening and discomfort. Bacteria from sweat and moisture can also build up over time.



Replacing your mask cushions on a regular basis is important for hygiene, which may help your therapy be as comfortable and effective as possible.



CPAP machine & accessories replacement tips

These replacement schedule guidelines are based on the current U.S. Medicare replacement schedule.*

☑ **Every month:** Disposable air filters

Filters can wear out with use. Dirt and dust can also build up over time, which can lead to blockages and unsanitary conditions.

☑ **Every 6 months:** Humidifier water chambers

Tubs can become discolored, cloudy, pitted or cracked, which can trap bacteria.

☑ **Every 3 months:** Air tubing

Tubing can develop small holes or tears over time, which can cause air leaks.

☑ **Every 5 years:** CPAP, APAP and Bilevel machine

You may need a new machine to enjoy new features that enhance your treatment and benefit from digital health tools offering personalized support for sleep therapy.

* The reimbursement information is being provided on an "as is" basis with no express or implied warranty of any kind and should be used solely for your internal informational purposes. Timelines provided are U.S. Medicare replacement timeframes. Co-pays, deductibles and co-insurance may apply. The information does not constitute professional or legal advice on reimbursement and should be used at your sole liability and discretion. All coding, coverage policies and reimbursement information are subject to change without notice. ResMed does not represent or warrant that any of the information being provided is true or correct and you agree to hold ResMed harmless in the event of any loss, damage, liabilities or claims arising from the use of the reimbursement information provided to you. Before filing any claims, it is the provider's sole responsibility to verify current requirements and coverage policies with the applicable payer.

Distributed by Resmed Corp, 9001 Spectrum Center Boulevard, San Diego, CA 92123 USA (858) 836-5000 or (800) 424-0737 (toll-free). See Resmed.com for other Resmed locations worldwide. Specifications may change without notice. © 2025 Resmed.