

## Style and design with her in mind



There's nothing hard about it. Swift™ FX Bella is so soft, light and simple, it feels like it's barely there! Minimal contact on the face and unobtrusive design make therapy less intimidating and easy to accept. It accommodates a variety of hairstyles too.

### Benefits for you ...

#### **Intuitive, almost instant fit**

Once the right pillows size is selected, Swift FX Bella requires minimal or often, no adjustment. So it saves you time, money and additional patient visits.

#### **Interchangeable with Swift FX for Her**

Swift FX Bella is fully interchangeable with Swift FX for Her (both use the same pillows and tubing) for an increased level of patient compliance.

#### **Easier therapy acceptance**

Designed for instant appeal, Swift FX Bella is not intimidating, helping users accept and continue therapy. And with very few parts, the system is easy to clean and assemble.

#### **Stable seal**

The Swift FX Bella system stays in place and maintains seal effectively.

### ... and your patients

#### **Hair management**

The exclusive Bella headgear is designed to stay out of her hair for uncomplicated hair management and the entire system enhances comfort.

#### **Interchangeable with Swift FX for Her**

It is fully interchangeable with Swift FX for Her (both have the same pillows and tubing) which provides users with another headgear option.

#### **Choice of sleeping positions**

Swift FX Bella is stable during movement as the user sleeps—so they can sleep on their side or their back, positioning the tube as desired.

#### **Quiet performance**

Innovative vent design diffuses airflow for quiet mask performance—only 25 dBA.\*

\* Declared dual numbered noise emission values (with uncertainty 3dBA) in accordance with ISO 4871.e 1



**Swift™ FX Bella**  
NASAL PILLOWS SYSTEM



## Caters for her personal style

Whether the user has long or short hair, Swift FX Bella is comfortable and convenient ... and it's finished in a stylish feminine colour scheme, just for her.

Pillows for her (XS, S and M) come packaged with the mask system to fit most female patients.

### Pillows for her...



**Flexible pillows with trampoline action**  
accommodate movement in multiple directions while maintaining seal

**Dual-wall nasal pillows**  
seal on contact for a gentle, comfortable seal, offer a wide fit range and disperse air gently to the nose

**Soft, narrow cushion with integrated flexible chamber**  
positions the tubing for ease of use

**Quiet vents**  
diffuse and direct airflow quietly away from the user and bed partner

## Product Codes

Swift FX Bella mask system *	61561
Swift FX Bella headgear only **	61581
Swift FX Bella headgear and Swift FX for Her headgear Combo **	61580

\* Available in Europe 1 only. Includes XS, S, and M pillows, Swift FX Bella headgear assembly, Swift FX for Her headgear assembly, pink backstrap and pink soft wraps. L pillows available separately.

\*\* Available in Europe 1, 2, 3.


**Europe 1:** English (UK), French, German and Italian (Switzerland)

**Europe 2:** English, Danish, Finnish, Norwegian and Swedish

**Europe 3:** English, Czech, Dutch (Belgium), French (Belgium), German (Belgium), Greek, Hungarian, Italian, Polish, Portuguese, Russian, Slovak, Spanish and Turkish



**Bella headgear**  
fits comfortably around the ears and stays out of her hair, for a more streamlined look and feel



**Ultra-smooth 360° rotating elbow**  
rotates with easy, fluid movement for flexibility in tube positions

**Spring-flex lightweight tube** stretches and bends uniquely for increased stability

**'Lock and key'**  
for easy attaching and detaching of headgear

**Size and name indicators**  
make it easy to identify the size and the correct orientation of the mask

**Swift™ FX Bella**  
NASAL PILLOWS SYSTEM