RESMED

Mirage Liberty™

FULL FACE MASK

Quick Fitting Guide



Components Card

English • Русский • Türkçe • Česky • Eesti • Slovensko • Български • Limba română • Slovensky • Lietuvių kalba • Latviski


Mirage Liberty Full Face Mask

Mask System / Система с маской / Maske sistemi / Systém masky / Maskikomplekt / Sistem maske / Устройство на маската / Sistem mască / Systém masky / Kaukės sistema / Maskas sistēma

	Europe 3*
S*	**61314**
L*	**61315**

Elbow Assembly (elbow, valve, clip) / Угловой патрубок в сборе (патрубок, клапан, зажим) / Dirsek takımı (dirsek, valf, klips) / Rameno sada (rameno, ventil, spojka) / Torupõlv (põlv, ventiil, klamber) / Kolenski sklop (koleno, ventil, objemka) / Комплект коляно (коляно, клапа, щипка) / Ansamblu cot (cot, valvă, clemă) / Kolenná jednotka (koleno, ventil, spona) / Alkūninės jungties mechanizmas (alkūninė jungtis, vožtuvas, vožtuvo spaustukas) / Izliekuma mezgls (izliekums, vārsts, fiksators)

61340



Europe 3*	English, Russian, Turkish, Czech, Estonian, Slovenian, Bulgarian, Romanian, Slovak, Lithuanian, Latvian
S*	Small / Маленький / Küçük / Malá / Väike / Majhna / Малък / Mică / Malá / Mažas / Mazs
M*	Medium / Средний / Orta / Střední / Keskmine / Srednja / Среден / Medie / Stredná / Vidutinis / Vidējs
L*	Large / Большой / Büyük / Velká / Suur / Velika / Голям / Mare / Veľká / Didelis / Liels
St*	Standard / Стандартный / Standart / Standardní / Standardne / Standardni / Стандартен / Standard / Štandardná / Standartiniai / Standarta

CE 0123

⚠ For details, refer to the *Mirage Liberty User Guide*. Specifications may change without notice. Mirage and Mirage Liberty are trademarks of ResMed Pty Ltd and Mirage is registered in U.S. Patent and Trademark Office. © 2011 ResMed. 618226/1 2011-11

618226